


Introducing

# ZB-DHA-PONA

- Designed and Tested for Detailed Hydrocarbon Analysis
- Excellent Response and Peak Symmetry for Polar Oxygenates
- Suitable for DHA, PONA, PIONA, PIANO, and ASTM Methods (D5134, D5441, D5501, D6729, D6730, D6733)

**Zebron™**  
GC Columns




**phenomenex®**  
...breaking with tradition™

[www.phenomenex.com/dhapona](http://www.phenomenex.com/dhapona)

# The Optimal Choice for PIONA, DHA and PONA

## Say Hello to the Zebron™ ZB-DHA-PONA GC Column

Zebron ZB-DHA-PONA is the optimal choice for the analysis of Detailed Hydrocarbon Analysis (DHA) within the Fuel industry. It's unique Engineered Self Cross-linking™(ESC) stationary phase provides low bleed, exceptional column life, and optimal separation of DHA critical pairs with symmetric peaks. In addition, Zebron ZB-DHA-PONA GC columns provide excellent response and peak symmetry for polar oxygenates.



### Easy ZB-DHA-PONA Selection for Your ASTM Method

Method	Description	Recommended Column	Recommended Dimensions	Phenomenex Part Number
ASTM D5134	Standard Test Method for Detailed Analysis of Petroleum Naphtha's through n-Nonane by Capillary Gas Chromatography	ZB-DHA-PONA	50 m x 0.20 mm x 0.5 µm	<a href="#">7JE-G042-17</a>
		ZB-DHA-PONA	100 m x 0.25 mm x 0.5 µm	<a href="#">7MG-G042-17</a>
ASTM D5441	Standard Test Method for Analysis of Methyl Tert-Butyl Ether (MTBE) by GC	ZB-DHA-PONA	50 m x 0.20 mm x 0.5 µm	<a href="#">7JE-G042-17</a>
		ZB-DHA-PONA	100 m x 0.25 mm x 0.5 µm	<a href="#">7MG-G042-17</a>
		ZB-DHA-PONA	150 m x 0.25 mm x 1 µm	<a href="#">7QG-G042-22</a>
ASTM D5501	Standard Test Method for Determination of Ethanol and Methanol Content in Fuels Containing Greater than 20% Ethanol by Gas Chromatography	ZB-DHA-PONA	100 m x 0.25 mm x 0.5 µm	<a href="#">7MG-G042-17</a>
		ZB-DHA-PONA	150 m x 0.25 mm x 1 µm	<a href="#">7QG-G042-22</a>
ASTM D6729	Standard Test Method for Determination of Individual Components in Spark Ignition Engine Fuels by 100 Meter Capillary High Resolution Gas Chromatography	ZB-DHA-PONA	100 m x 0.25 mm x 0.5 µm	<a href="#">7MG-G042-17</a>
ASTM D6730	Standard Test Method for Determination of Individual Components in Spark Ignition Engine Fuels by 100-Meter Capillary (with Pre-column) High-Resolution Gas Chromatography	ZB-DHA-PONA	50 m x 0.20 mm x 0.5 µm	<a href="#">7JE-G042-17</a>
		ZB-DHA-PONA	100 m x 0.25 mm x 0.5 µm	<a href="#">7MG-G042-17</a>
		ZB-DHA-PONA	150 m x 0.25 mm x 1 µm	<a href="#">7QG-G042-22</a>
ASTM D6733	Standard Test Method for Determination of Individual Components in Spark Ignition Engine Fuels by 50-Meter Capillary High Resolution Gas Chromatography	ZB-DHA-PONA-TUNE	5 m x 0.25 mm x 1 µm	<a href="#">7AG-G042-22</a>
		ZB-DHA-PONA	50 m x 0.20 mm x 0.5 µm	<a href="#">7JE-G042-17</a>

## Discover the Zebron™ ZB-DHA-PONA

Detailed Hydrocarbon Mixture Analysis: PIONA .....	4-5
Analysis of MTBE and Related Compounds: ASTM D5441 .....	6
Improved Peak Shape and Resolution for ASTM D5441.....	7
Determination of Individual Components in Spark Ignition Engine Fuels: ASTM D6730 .....	8
Quality Tested to Ensure Exceptional Column Reproducibility .....	9

## Additional GC Applications for Fuel Analysis

Simulated Distillation: Zebron ZB-1XT SimDist .....	10
High Temperature Hydrocarbon Analysis: Zebron ZB-HT Inferno.....	11

## Further Support for Your Analysis

Zebron Gas Management Filters.....	12
Zebron Gas Management Click-on Traps.....	13
Ordering Information.....	14-15



# PIONA Analysis

We developed this special Engineered Self Cross-linking™ (ESC) stationary phase with an optimized 100% dimethylpolysiloxane to allow for an ideal separation and resolution of all the components needed for PIONA and Detailed Hydrocarbon Analysis (DHA). The ESC technology found in each Zebron ZB-DHA-PONA also provides the additional benefit of low bleed and exceptional column life.

- Highly efficient dimension (100 meter length x 0.25 mm ID x 0.50 μm film thickness) provides ideal separation and resolution
- The smaller 0.25 mm diameter along with the column length provides an efficient mass transfer of the complex analytes for great separation and peak shapes
- 100% dimethylpolysiloxane provides true boiling point-based separation

The PIONA analysis below illustrates the successful separation of Paraffins, Isoparaffins, Olefins, Naphthenes, and Aromatics using one ZB-DHA-PONA GC column dimension and one single method.


### Conditions for all applications:

**Column:** Zebron™ ZB-DHA-PONA  
**Dimensions:** 100 meter x 0.25 mm x 0.50 μm  
**Part No.:** ZMG-G042-17  
**Injection:** Split 40:1 @ 300 °C, 0.2 μL

**Recommended Liner:** Zebron PLUS Straight Z-Liner™

**Part No.:** AG2-0A03-05  
**Carrier Gas:** Helium @ 1.55 mL/min (Constant Flow)  
**Oven Program:** 35 °C for 14 min, to 60 °C @ 1.1 °C/min for 19 min, to 280 °C @ 2 °C/min for 5 min  
**Detector:** FID @ 320 °C


## Separation of Paraffins



### Sample:

1. Pentane
2. Hexane
3. Heptane
4. Octane
5. Nonane
6. Decane
7. Undecane
8. Dodecane
9. Tridecane
10. Tetradecane
11. Pentadecane


## Separation of Isoparaffins



### Sample:

- |                        |                          |                            |                         |                         |                        |
|------------------------|--------------------------|----------------------------|-------------------------|-------------------------|------------------------|
| 1. Isopentane          | 7. 2,2,3-Trimethylbutane | 13. 2,2-Dimethyl hexane    | 19. 4-Methylheptane     | 25. 3,5-Dimethylheptane | 31. 3,3-Dimethyloctane |
| 2. 2,3-Dimethylbutane  | 8. 3,3-Dimethylpentane   | 14. 2,5-Dimethyl hexane    | 20. 3-Methylheptane     | 26. 3,4-Dimethylheptane | 32. 2,3-Dimethyloctane |
| 3. 2-Methylpentane     | 9. 2-Methylhexane        | 15. 2,2,3-Trimethylpentane | 21. 3-Ethylhexane       | 27. 2-Methyloctane      | 33. 3-Ethylhexane      |
| 4. 3-Methylpentane     | 10. 2,3-Dimethylpentane  | 16. 2,4-Dimethylhexane     | 22. 2,5-Dimethylheptane | 28. 3,3-Diethylpentane  | 34. 2-Methylnonane     |
| 5. 2,2-Dimethylpentane | 11. 3-Methylhexane       | 17. 2,3-Dimethylhexane     | 23. 3,3-Dimethylheptane | 29. 3-Methyloctane      | 35. 3-Methylnonane     |
| 6. 2,4-Dimethylpentane | 12. 3-Ethylpentane       | 18. 2-Methylheptane        | 24. 2,3-Dimethylheptane | 30. 2,2-Dimethyloctane  |                        |


## Separation of Aromatic Compounds



**Sample:**

- |                     |                             |                                 |                                   |                                   |
|---------------------|-----------------------------|---------------------------------|-----------------------------------|-----------------------------------|
| 1. Benzene          | 9. 1-Methyl-3-ethylbenzene  | 17. 1-Methyl-3-isopropylbenzene | 25. 1,4-Dimethyl-2-ethylbenzene   | 33. n-Pentylbenzene               |
| 2. Toluene          | 10. 1-Methyl-4-ethylbenzene | 18. 1-Methyl-4-isopropylbenzene | 26. 1,3-Dimethyl-5-ethylbenzene   | 34. t-1-Butyl,3,5-dimethylbenzene |
| 3. Ethylbenzene     | 11. 1,3,5-Trimethylbenzene  | 19. 1-Methyl-2-isopropylbenzene | 27. 1,2-Dimethyl-4-ethylbenzene   | 35. t-1-butyl-ethylbenzene        |
| 4. m-Xylene         | 12. 1-Methyl-2-ethylbenzene | 20. 1-Methyl-3-n-propylbenzene  | 28. 1,3-Dimethyl-2-ethylbenzene   | 36. 1,3,5-Triethylbenzene         |
| 5. p-Xylene         | 13. 1,2,4-Trimethylbenzene  | 21. 1-Methyl-4-n-propylbenzene  | 29. 1,2-Dimethyl-3-ethylbenzene   | 37. 1,2,4-Triethylbenzene         |
| 6. o-Xylene         | 14. tert-Butylbenzene       | 22. n-Butylbenzene              | 30. 1,2,4,5-Tetramethylbenzene    | 38. n-Hexylbenzene                |
| 7. Isopropylbenzene | 15. Isobutylbenzene         | 23. 1-Methyl-2-n-propylbenzene  | 31. 2-Methylbutylbenzene          |                                   |
| 8. n-Propylbenzene  | 16. sec-Butylbenzene        | 24. 1,2-Diethylbenzene          | 32. trans-1-Butyl-2-methylbenzene |                                   |


## Separation of Naphthenes



**Sample:**

- |                                   |                                     |                                     |                                    |                                    |
|-----------------------------------|-------------------------------------|-------------------------------------|------------------------------------|------------------------------------|
| 1. Cyclopentane                   | 7. trans-1,2-Dimethylcyclopentane   | 13. trans-1,4-Dimethylcyclohexane   | 19. n-Propylcyclopentane           | 25. Isobutylcyclopentane           |
| 2. Methylcyclopentane             | 8. Methylcyclohexane                | 14. 1-Ethyl-1-methylcyclopentane    | 20. ccc-1,3,5-Trimethylcyclohexane | 26. Isopropylcyclohexane           |
| 3. Cyclohexane                    | 9. Ethylcyclopentane                | 15. trans-1,2-Dimethylcyclohexane   | 21. 1,1,4-Trimethylcyclohexane     | 27. n-Butylcyclopentane            |
| 4. 1,1-Dimethylcyclopentane       | 10. ctc-1,2,4-Trimethylcyclopentane | 16. ccc-1,2,3-Trimethylcyclopentane | 22. ctt-1,2,4-Trimethylcyclohexane | 28. Isobutylcyclohexane            |
| 5. cis-1,3-Dimethylcyclopentane   | 11. ctc-1,2,3-Trimethylcyclopentane | 17. Isopropylcyclopentane           | 23. ctc-1,2,4-Trimethylcyclohexane | 29. t-1-Methyl-2-propylcyclohexane |
| 6. trans-1,3-Dimethylcyclopentane | 12. cct-1,2,4-Trimethylcyclopentane | 18. cis-1,2-Dimethylcyclohexane     | 24. 1,1,2-Trimethylcyclohexane     | 30. t-1-Methyl-2-(4MP)cyclopentane |

## Separation of Olefins



**Sample:**

- |                           |                      |                     |                    |                    |
|---------------------------|----------------------|---------------------|--------------------|--------------------|
| 1. 3-Methyl-1-butene      | 6. cis-2-Pentene     | 11. cis-2-Hexene    | 16. cis-2-Heptene  | 21. trans-3-Nonene |
| 2. 1-Pentene              | 7. 4-Methylpentene-1 | 12. 1-Heptene       | 17. 1-Octene       | 22. cis-3-Nonene   |
| 3. 2-Methyl-1-butene      | 8. 1-Hexene          | 13. trans-3-Heptene | 18. trans-2-Octene | 23. trans-2-Nonene |
| 4. 2-Methyl-1,3-butadiene | 9. 2-Methylpentene-2 | 14. cis-3-Heptene   | 19. cis-2-Octene   | 24. cis-2-Nonene   |
| 5. trans-2-Pentene        | 10. trans-2-Hexene   | 15. trans-2-Heptene | 20. 1-Nonene       | 25. 1-Decene       |

# ASTM D5441


## Analysis of MTBE and Related Impurities: ASTM D5441

ASTM D5441 is a gas chromatography analysis used to measure Methyl-t-Butyl Ether (MTBE) as its purity is critical for the performance of spark ignition fuel.

Zebtron ZB-DHA-PONA easily identifies impurities left during the manufacturing process of MTBE and that no contaminants were introduced in its storage or shipment. This ASTM D5441 analysis demonstrates both optimal peak shape and resolution with a Zebtron ZB-DHA-PONA GC column:

- Achieve symmetric peaks for oxygenates like MTBE
- Meet and exceed ASTM D5441 resolution requirements

### Analysis of MTBE and related compounds on Zebtron ZB-DHA-PONA by GC-FID



### Zebtron ZB-DHA-PONA Exceeds ASTM D5441 Resolution Requirements


Analyte	Resolution	
	ASTM D5441 Requirement	Zebtron ZB-DHA-PONA
Trans-2-Pentene and Tert-Butanol	Not less than 1.3	1.6
Tert-Butanol and Cis-2-Pentene	Not less than 1.3	1.5

## Improved Peak Shape and Resolution for ASTM D5441


By switching to a Zebron ZB-DHA-PONA for the analysis of Methyl-t-Butyl Ether (MTBE) and related compounds, you gain optimal separation of DHA critical pairs with symmetric peak shape.

Analyte Peaks	ASTM Requirement	Resolution			
		Zebron ZB-DHA-PONA	Restek® Rtx®-DHA	Supelco® Petrocol® DH 50.2	Agilent® HP®-PONA
Trans-2-Pentene and Tert-Butanol	Not less than 1.3	1.6	1.4	2.1	1.7
Tert-Butanol and Cis-2-Pentene	Not less than 1.3	1.5	1.5	0.8	1.2


### Zebron ZB-DHA-PONA




### Restek Rtx-DHA



### Supelco Petrocol DH 50.2



### Agilent HP-PONA



#### GC Conditions for all columns:

- Columns: As noted
- Dimensions: 50 meter x 0.20 mm x 0.50  $\mu$ m
- Part No.: [7JE-G042-17](#)
- Injection: Pulse Split 50:1, 20 psi for 0.5 min, @ 200 °C, 0.5  $\mu$ L
- Recommended Liner: Zebron PLUS Z-Liner™
- Part No.: [AG2-0A03-05](#)
- Carrier Gas: Helium @ 0.9 mL/min (Constant Flow)
- Oven Program: 40 °C for 13 min, to 180 °C at 10 °C/min for 3 min
- Detector: FID @ 250 °C

Comparative separations may not be representative of all applications.


# ASTM D6730

## Determination of Individual Components in Spark Ignition Engine Fuels

ASTM D6730 is an important standard test method for the determination of individual components in Spark Ignition Engine Fuels. By using a 100 meter Zebtron ZB-DHA-PONA GC column in conjunction with the ZB-DHA-PONA-TUNE tuning column, you can achieve High-Resolution Gas Chromatography separation of necessary critical pairs.


**Column 1 (Tuning):** Zebtron ZB-DHA-PONA-TUNE  
**Dimensions:** 5 meter x 0.25 mm x 1.0 μm  
**Part No.:** ZAG-G042-22  
**Column 2:** Zebtron ZB-DHA-PONA  
**Dimensions:** 100 meter x 0.25 mm x 0.50 μm  
**Part No.:** ZMG-G042-17

**Recommended Column Union:** AG0-4716  
**Injection:** Split 150:1 @ 200 °C, 0.2 μL  
**Recommended Liner:** Zebtron PLUS Straight Z-Liner™  
**Part No.:** AG2-0A03-05  
**Carrier Gas:** Hydrogen @ 2 mL/min (Constant Flow)  
**Oven Program:** 30 °C for 8.5 min, to 48 °C @ 22 °C/min. for 27 min, to 141 °C @ 3 °C/min for 1 min, to 275 °C @ 1 °C/min for 2 min  
**Detection:** FID @ 275 °C



### Tips

- Use a ZB-DHA-PONA-TUNE with ZB-DHA-PONA to get good separation of aromatics like benzene and toluene from paraffins.
- Additional ZB-DHA-PONA dimensions such as 50m and 150m are available for faster or higher efficiency analyses.




[www.phenomenex.com/dha-pona](http://www.phenomenex.com/dha-pona)



# Exceptional Column Reproducibility

## QC Tested for the Compounds You Analyze

We added challenging and troublesome analytes to our QC test to make sure each Zebtron ZB-DHA-PONA column provides superior deactivation, improved peak shape, and exhibits high efficiency and exceptional inertness, for nonpolar hydrocarbons and polar oxygenates in gasoline, so you can be sure your column is ready to meet suitability requirements for your methods.



Column	Retention Factor (k)						
	Undecane	4-Chlorophenol	Tridecane	1-Undecanol	Tetradecane	Dicyclohexylamine	Pentadecane
1	1.11	1.5	2.91	3.80	4.69	5.12	7.56
2	1.11	1.51	2.92	3.82	4.71	5.14	7.58
3	1.10	1.49	2.89	3.77	4.66	5.08	7.5
% RSD	0.52	0.67	0.53	0.66	0.54	0.60	0.55

Low % RSD indicates high reproducibility

Column	Peak Skew						
	Undecane	4-Chlorophenol	Tridecane	1-Undecanol	Tetradecane	Dicyclohexylamine	Pentadecane
1	1.0	1.0	1.0	1.0	1.0	1.0	1.0
2	1.0	1.0	1.0	1.0	0.9	0.9	0.9
3	1.0	1.0	1.0	1.0	0.9	0.9	1

Acceptable range for skew is 0.80 to 1.5. Zebtron ZB-DHA-PONA provides symmetric peaks for polar compounds (highlighted in blue) and nonpolar compounds (highlighted in yellow) which is important for DHA analysis.


# Zebtron ZB-1XT SimDist GC Column

## Simulated Distillation


Zebtron ZB-1XT SimDist GC Columns are designed to meet the rigorous requirements of ASTM, EN, and other specialized methods for the fuels industry.

- Uniform Glass Infusion™ Technology Coating provides sharper peaks and higher efficiency.
- Columns are individually tested to ensure excellent reproducibility and overall performance for SimDist critical compounds.
- 45-70% Higher column efficiency than traditional metal GC columns.
- Improved lifetime and resolution of C50/C52 hour after hour.
- 100% Risk-Free Guarantee.

### Glass Infused ZB-1XT



### Not Glass Infused



**Conditions for both columns:**

**Injection:** Split 100:1 @ 250 °C, 1.0 µL  
**Carrier Gas:** Hydrogen @ 1.93 mL/min (constant pressure)  
**Oven Program:** 130 °C (Isothermal)  
**Detector:** FID @ 325 °C

**Sample:** 1. Undecane  
 2. 4-Chlorophenol  
 3. Tridecane  
 4. 1-Undecanol  
 5. Dicyclohexylamine  
 6. Pentadecane

Note: Sample at 250 µg/mL each in hexane

## Zebtron ZB-1XT SimDist Column Selection per ASTM method

ASTM Method	Ranges	Applications	Recommended Column	Description
D7169	<a href="#">C9-C100</a>	Crude Oil	<a href="#">7AK-G026-05</a>	ZB-1XT SimDist, 5 m x 0.53 mm x 0.15 µm
D7169	<a href="#">C9-C100</a>	Vacuum Distillates	<a href="#">7AK-G026-55</a>	ZB-1XT SimDist, 5 m x 0.53 mm x 0.09 µm
D6352	C10-C90	Petroleum Distillates	<a href="#">7AK-G026-05</a>	ZB-1XT SimDist, 5 m x 0.53 mm x 0.15 µm
D2887X	C10-C60	Petroleum Fractions	<a href="#">7AK-G026-05</a>	ZB-1XT SimDist, 5 m x 0.53 mm x 0.15 µm
D2887	C5-C44	Petroleum Fractions	<a href="#">7CK-G026-35</a>	ZB-1XT SimDist, 10 m x 0.53 mm x 2.65 µm
D3710	C4-C15	Gasoline Fractions	<a href="#">7CK-G026-35</a>	ZB-1XT SimDist, 10 m x 0.53 mm x 2.65 µm

# Zebron ZB-1HT Inferno™ GC Column

## High Temperature Hydrocarbon Analysis


### Rugged GC Performance to 430°C!

Zebron Inferno columns allow you to push your standard analysis up to 430°C for increased lab productivity and improved performance – especially if you struggle with high boiling compounds, contaminants, or carry-overs.

Zebron Inferno columns have the inertness and the temperature stability needed to separate high molecular weight compounds. These high boilers require a thermally rugged column with low bleed and able to provide optimal resolution like the ZB-1HT.

- Better peak shape
- Improved separation
- More analytes detected

### Great Separation of High Boiling Hydrocarbons (ASTM Method D6352)



#### Easily Elute Hydrocarbons Up To C90

Zebron ZB-1HT Inferno can operate up to 430°C, allowing C90 to elute on temperature ramp to meet ASTM requirements!

#### Sample:

1. C10	13. C20	25. C32	37. C48	49. C72
2. C11	14. C21	26. C33	38. C50	50. C74
3. C12	15. C22	27. C34	39. C52	51. C76
4. C13	16. C23	28. C35	40. C54	52. C78
5. C14	17. C24	29. C36	41. C56	53. C80
6. C15	18. C25	30. C37	42. C58	54. C82
7. C16	19. C26	31. C38	43. C60	55. C84
8. C17	20. C27	32. C39	44. C62	56. C86
9. Pristane	21. C28	33. C40	45. C64	57. C88
10. C18	22. C29	34. C42	46. C66	58. C90
11. Phytane	23. C30	35. C44	47. C68	
12. C19	24. C31	36. C46	48. C70	

Note: Sample was a combination of POLYWAX® 655 and retention time markers C8-C40 in CS<sub>2</sub>/Chloroform

ZB-1HT Inferno

**Column:** Zebron ZB-1HT Inferno  
**Dimensions:** 5 meter x 0.53 mm x 0.10 µm  
**Part No.:** [ZAK-G014-02](#)  
**Injection:** On-Column @ 43°C, 0.1 µL  
**Carrier Gas:** Helium @ 4.4 mL/min (constant flow)  
**Oven Program:** 40°C for 0.5 min to 430°C @ 20°C/min for 10 min  
**Detector:** FID @ 430°C



# Zebtron Gas Management Filters

## For Ultra Pure GC Gas Supply

### No Diffusion and No Leaks with Dual Wall Housing!

More protection and less contamination with an inner glass wall and better safety with a plastic outer wall.

### Higher Capacity

The higher capacity of Zebtron gas management filters means less need for replacement and an increase in system uptime, saving you time and money.

### Easy Filter Installation in Under 20 Seconds!

Avoid system downtime by attaching a new filter to your gas line in seconds!

### Never Disrupt Your Instrument!

The integrated one way valve on the Zebtron gas management base units allow for quick and easy filter changes without disrupting the Instrument.



### Easily Know When to Replace Your Filter!

1

#### Color Change

The indicator on each Zebtron gas management filter will display when it is time to change your filter.\*

2

#### Electronic Indicator

The optional electronic indicator will alert you with an audible warning on when to change your filter based on typical usage.

### Universal Compatibility

No need to change your existing connecting unit. Simply remove the old filter, replace the O-ring, click a new filter in, and tighten the connection with the ring nut! (Part No.: [AG6-1060](#)).

\*Hydrocarbon filters do not include a color indicator

# Zebron Gas Management Traps

## For Ultra Pure GC Gas Supply

### Click-on Trap

Never Disrupt Your Instrument!  
Attach a new filter in under 10 seconds. Once the connecting unit is installed, you will be able to change traps as often as needed without interrupting the instrument.

### High Pressure Durability

The Click-on trap is completely made out of stainless steel and can withstand high pressures. The seals at each end of the trap will only be punctured when the Click-on connector is locked in place.

### Know When to Replace

The optional electronic indicator displays when scheduled replacement or maintenance is due.

### Secure Wall Mount

These wall mounting clamps will EASILY help attach a Zebron gas trap to a wall or surface.



### Click-On Connection

The Click-on trap connectors allow inline cartridges to be exchanged without introducing contaminants. Spring loaded check valves seal when a filter is removed and open only when a new filter has been locked in place.



# Gas Management Ordering Information



## Zebtron Gas Management Filters

Part No.	Description	Unit	Price
<a href="#">AG6-1040</a>	Gas Filter Universal	ea	\$ 185
<a href="#">AG6-1070</a>	Gas Filter Universal (Helium)	ea	210
<a href="#">AG6-1010</a>	Gas Filter Oxygen	ea	160
<a href="#">AG6-1020</a>	Gas Filter Moisture	ea	160
<a href="#">AG6-1030</a>	Gas Filter Hydrocarbon	ea	160
<a href="#">AG6-1050</a>	Gas Filter Hydrocarbon/Moisture for LC-MS	2/pk	211
<a href="#">AG6-1060</a>	Ring Nut for Gas Filter	ea	16



## Zebtron Gas Management Traps

Part No.	Description	Unit	Price
<a href="#">AG6-3140</a>	Click-on Universal Trap	ea	\$ 338
<a href="#">AG6-3110</a>	Click-on Oxygen Trap	ea	256
<a href="#">AG6-3120</a>	Click-on Moisture Trap	ea	294
<a href="#">AG6-3130</a>	Click-on Hydrocarbon Trap	ea	294
<a href="#">AG6-3150</a>	Click-on Carbon Dioxide Trap	ea	307



## Zebtron Connecting Units

Part No.	Description	Unit	Price
<a href="#">AG6-2101</a>	1-position Connecting Unit 1/4 in. Brass	ea	\$ 206
<a href="#">AG6-2102</a>	2-position Connecting Unit 1/4 in. Brass	ea	391
<a href="#">AG6-2103</a>	4-position Connecting Unit 1/4 in. Brass	ea	747
<a href="#">AG6-2201</a>	1-position Connecting Unit 1/8 in. Brass	ea	206
<a href="#">AG6-2202</a>	2-position Connecting Unit 1/8 in. Brass	ea	391
<a href="#">AG6-2203</a>	4-position Connecting Unit 1/8 in. Brass	ea	747
<a href="#">AG6-2204</a>	High flow 2-position Connecting Unit for LC-MS	ea	697
<a href="#">AG6-2205</a>	Particle Filter for LC-MS	ea	98
<a href="#">AG6-2206</a>	O-ring Replacement for Gas Filter Baseplate	20/pk	29

Part No.	Description	Unit	Price
<a href="#">AG6-2301</a>	1-position Connecting Unit 1/4 in. Stainless Steel	ea	\$ 320
<a href="#">AG6-2302</a>	2-position Connecting Unit 1/4 in. Stainless Steel	ea	600
<a href="#">AG6-2303</a>	4-position Connecting Unit 1/4 in. Stainless Steel	ea	950
<a href="#">AG6-2304</a>	1-position Connecting Unit 1/8 in. Stainless Steel	ea	320
<a href="#">AG6-2305</a>	2-position Connecting Unit 1/8 in. Stainless Steel	ea	600
<a href="#">AG6-2306</a>	4-position Connecting Unit 1/8 in. Stainless Steel	ea	950

## Zebtron Base Electronic Indicator and Other Accessories

Part No.	Description	Unit	Price
<a href="#">AG6-3160</a>	1/8 in. Brass Click-on Connector Set	2/pk	\$ 39
<a href="#">AG6-3170</a>	1/4 in. Brass Click-on Connector Set	2/pk	56
<a href="#">AG6-3180</a>	Wall-mounting Clamp Set for Gas Traps	2/pk	150
<a href="#">AG6-4150</a>	1/8 in. Stainless Steel Click-on Connector Set	2/pk	60
<a href="#">AG6-4160</a>	1/4 in. Stainless Steel Click-on Connector Set	2/pk	70
<a href="#">AG6-3190</a>	O-ring Replacement Set for Gas Trap	20/pk	25
<a href="#">AG6-4110</a>	Electronic Indicator for Gas Trap	ea	60
<a href="#">AG6-4120</a>	Electronic Indicator for Gas Filter	ea	60
<a href="#">AG6-4130</a>	Electronic Indicator for LC-MS Filter	ea	60



## First time ordering Zebtron Click-on Gas Traps?

Be sure to order the **brass or stainless steel connector** with your first trap.



# GC Column

## Ordering Information


### Zebron ZB-DHA-PONA

ID (mm)	df (µm)	Temp. Limits °C	Part No.	Price
<b>50-Meter</b>				
0.2	0.5	-60 to 360/370C	<a href="#">7JE-G042-17</a>	\$ 799
<b>100-Meter</b>				
0.25	0.5	-60 to 360/370C	<a href="#">7MG-G042-17</a>	\$ 950
<b>150-Meter</b>				
0.25	1	-60 to 340/360C	<a href="#">7QG-G042-22</a>	\$ 1,400



### Zebron ZB-DHA-PONA-TUNE

ID (mm)	df (µm)	Temp. Limits °C	Part No.	Price
<b>5-Meter</b>				
0.25	1	-60 to 340/360C	<a href="#">7AG-G042-22</a>	\$ 95



### Zebron ZB-1XT SimDist

ID(mm)	df (µm)	Temp. Limits °C	Part No.	Price
<b>5-Meter</b>				
0.53	0.09	-60 to 450	<a href="#">7AK-G026-55</a>	\$ 430
0.53	0.15	-60 to 450	<a href="#">7AK-G026-05</a>	430
<b>5-Meter with 2-Meter Guardian™ Integrated Guard</b>				
0.53	0.09	-60 to 450	<a href="#">7AK-G026-55-GGT</a>	\$ 450
0.53	0.15	-60 to 450	<a href="#">7AK-G026-05-GGT</a>	450
<b>10-Meter</b>				
0.53	0.15	-60 to 450	<a href="#">7CK-G026-05</a>	\$ 425
0.53	0.88	-60 to 450	<a href="#">7CK-G026-49</a>	600
0.53	2.65	-60 to 450	<a href="#">7CK-G026-35</a>	490
<b>15-Meter</b>				
0.53	0.25	-60 to 450	<a href="#">7EK-G026-11</a>	\$ 530



If you need a 5 in. cage, simply add a (-B) after the part number, e.g., [7CK-G026-05-B](#). Some exceptions may apply. Agilent 6850 and some SRI and process GC systems use only 5 in. cages.

### Zebron ZB-1HT Inferno

ID (mm)	df (µm)	Temp. Limits °C	Part No.	Price
<b>5-Meter</b>				
0.53	0.10	-60 to 400/430C	<a href="#">7AK-G014-02</a>	\$ 350



To Order Zebron GC Columns, go to [www.phenomenex.com/Zebron](http://www.phenomenex.com/Zebron)

**BE-HAPPY™**  
guarantee

Your happiness is our mission. Take 45 days to try our products. If you are not happy, we'll make it right.

[www.phenomenex.com/behappy](http://www.phenomenex.com/behappy)

## Zebron PLUS Liners

We've made it easier than ever to find your ideal GC liner. Search by:

- Application
- GC System
- Injection Type
- Your Current Liner Part Number

[www.phenomenex.com/findLiner](http://www.phenomenex.com/findLiner)



# Zebron ZB-DHA-PONA

The Detailed Hydrocarbon Analysis (DHA) GC column designed specifically for DHA, PONA, PIONA, and PIANO analysis

**Australia**

t: +61 (0)2-9428-6444  
auiinfo@phenomenex.com

**Austria**

t: +43 (0)1-319-1301  
anfrage@phenomenex.com

**Belgium**

t: +32 (0)2 503 4015 (French)  
t: +32 (0)2 511 8666 (Dutch)  
beinfo@phenomenex.com

**Canada**

t: +1 (800) 543-3681  
info@phenomenex.com

**China**

t: +86 400-606-8099  
cninfo@phenomenex.com

**Denmark**

t: +45 4824 8048  
nordicinfo@phenomenex.com

**Finland**

t: +358 (0)9 4789 0063  
nordicinfo@phenomenex.com

**France**

t: +33 (0)1 30 09 21 10  
franceinfo@phenomenex.com

**Germany**

t: +49 (0)6021-58830-0  
anfrage@phenomenex.com

**India**

t: +91 (0)40-3012 2400  
indiainfo@phenomenex.com

**Ireland**

t: +353 (0)1 247 5405  
eireinfo@phenomenex.com

**Italy**

t: +39 051 6327511  
italiainfo@phenomenex.com

**Luxembourg**

t: +31 (0)30-2418700  
nlinfo@phenomenex.com

**Mexico**

t: 01-800-844-5226  
tecnicomx@phenomenex.com

**The Netherlands**

t: +31 (0)30-2418700  
nlinfo@phenomenex.com

**New Zealand**

t: +64 (0)9-4780951  
nzinfo@phenomenex.com

**Norway**

t: +47 810 02 005  
nordicinfo@phenomenex.com

**Portugal**

t: +351 221 450 488  
ptinfo@phenomenex.com

**Singapore**

t: +65 800-852-3944  
sginfo@phenomenex.com

**Spain**

t: +34 91-413-8613  
espinfo@phenomenex.com

**Sweden**

t: +46 (0)8 611 6950  
nordicinfo@phenomenex.com

**Switzerland**

t: +41 (0)61 692 20 20  
swissinfo@phenomenex.com

**Taiwan**

t: +886 (0) 0801-49-1246  
twinfo@phenomenex.com

**United Kingdom**

t: +44 (0)1625-501367  
ukinfo@phenomenex.com

**USA**

t: +1 (310) 212-0555  
info@phenomenex.com

**All other countries  
Corporate Office USA** 

t: +1 (310) 212-0555  
info@phenomenex.com



[www.phenomenex.com](http://www.phenomenex.com)

Phenomenex products are available worldwide. For the distributor in your country, contact Phenomenex USA, International Department at [international@phenomenex.com](mailto:international@phenomenex.com)

**Terms and Conditions**

Subject to Phenomenex Standard Terms and Conditions, which may be viewed at [www.phenomenex.com/TermsAndConditions](http://www.phenomenex.com/TermsAndConditions).

**Trademarks**

Engineered Self Cross-linking (ESC), Inferno, Zebron, and Z-Liner and are trademarks of Phenomenex. Agilent and HP are registered trademarks of Agilent Technologies, Inc. Restek and Rtx are registered trademarks of Restek Corporation. Supelco and Petrocol are registered trademarks of Sigma-Aldrich Co. LLC.

**Disclaimers**

Comparative separations may not be representative of all applications. Phenomenex is not affiliated with Agilent Technologies, Restek Corporation, or Sigma-Aldrich Co..

FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures.

© 2019 Phenomenex, Inc. All rights reserved.