

Tips from our Protein Separation ZenMasters

Drug-Antibody Ratio (DAR)



Chad Eichman, Ph.D.
BioPharm Global Marketing Manager


Why do average DAR for cysteine conjugates use Reversed Phase HPLC?

Average DAR is a critical quality attribute for Antibody-Drug Conjugate (ADC), as it gives indications of toxicity, efficacy, and pharmacokinetics.¹ As such, DAR is determined very early, even during early development research.


Although Hydrophobic Interaction Chromatography (HIC) might give both average DAR and drug distribution, this approach requires rather extensive mobile phase optimization and method development, and might be more appropriate further downstream of drug development. Another option to assess DAR for cysteine conjugates is Native MS. Although this approach gives assessment of drug distribution and DAR, it requires a high-resolution instrument such as a Time-of-Flight (TOF)-MS, not to mention expertise in large molecule MS.

Given the above challenges, reversed phase HPLC is a preferred route for average DAR. This approach can be done with cysteine conjugates simply by fully reducing the protein and summing together the total peak areas for conjugated heavy and light chains, as indicated below. With optimized conditions, there is also partial separation of different positional isomers, though again, drug distribution would need to be assessed with an orthogonal method.

Average DAR (Trastuzumab mc-PAB-MMAF Cysteine Conjugate)



Average DAR (Trastuzumab vc-PAB-MMAE Cysteine Conjugate)



1. Peters, C.; Brown, S. "Antibody-Drug conjugates as novel anti-Cancer chemotherapeutics." Bioscience Reports, vol. 35, no. 4, Dec. 2015, doi:10.1042/bsr20150089.

Find the conditions online at:
www.phenomenex.com/bioZen

Product Ordering Information

bioZen™ Products - Powered by Biocompatible Hardware

bioZen Columns (mm)	Phases							Biocompatible Guard Cartridges*	
	50 x 2.1	100 x 2.1	150 x 2.1	250 x 2.1	50 x 4.6	150 x 4.6	300 x 4.6	for 2.1 mm	for 4.6 mm
bioZen 2.6 µm Glycan	—	00D-4773-AN	00F-4773-AN	—	—	—	—	AJO-9800	—
bioZen 1.6 µm Peptide PS-C18	00B-4770-AN	—	00F-4770-AN	—	—	—	—	AJO-9803	—
bioZen 3 µm Peptide PS-C18	—	—	—	—	00B-4771-E0	00F-4771-E0	—	—	AJO-7606
bioZen 1.7 µm Peptide XB-C18	00B-4774-AN	—	00F-4774-AN	—	—	—	—	AJO-9806	—
bioZen 2.6 µm Peptide XB-C18	00B-4768-AN	—	00F-4768-AN	00G-4768-AN	00B-4768-E0	00F-4768-E0	—	AJO-9806	AJO-9808
bioZen 3.6 µm Intact C4	00B-4767-AN	—	00F-4767-AN	—	00B-4767-E0	00F-4767-E0	—	AJO-9809	AJO-9811
bioZen 3.6 µm Intact XB-C8	00B-4766-AN	—	00F-4766-AN	—	00B-4766-E0	00F-4766-E0	—	AJO-9812	AJO-9814
bioZen 1.8 µm SEC-2	—	—	—	—	—	00F-4769-E0	00H-4769-E0	—	AJO-9850
bioZen 1.8 µm SEC-3	—	—	—	—	—	00F-4772-E0	00H-4772-E0	—	AJO-9851

*AJO-7606 requires guard holder KJO-4282. All other guard cartridges require guard holder AJO-9000.



Have Questions? Visit:
www.phenomenex.com/bioZen
 We're here to help!

guarantee

If bioZen columns in this technical note do not provide at least equivalent separations as compared to a competing column of the same phase, particle size, and dimensions, return the column with the comparative data within 45 days for a FULL REFUND.

Australia

t: +61 (0)2-9428-6444
 f: +61 (0)2-9428-6445
 auiinfo@phenomenex.com

Austria

t: +43 (0)1-319-1301
 f: +43 (0)1-319-1300
 anfrage@phenomenex.com

Belgium

t: +32 (0)2 503 4015 (French)
 t: +32 (0)2 511 8666 (Dutch)
 f: +31 (0)30-2383749
 beinfo@phenomenex.com

Canada

t: +1 (800) 543-3681
 f: +1 (310) 328-7768
 info@phenomenex.com

China

t: +86 400-606-8099
 f: +86 (0)22 2532-1033
 phen@agela.com

Denmark

t: +45 4824 8048
 f: +45 4810 6265
 nordicinfo@phenomenex.com

Finland

t: +358 (0)9 4789 0063
 f: +45 4810 6265
 nordicinfo@phenomenex.com

France

t: +33 (0)1 30 09 21 10
 f: +33 (0)1 30 09 21 11
 franceinfo@phenomenex.com

Germany

t: +49 (0)6021-58830-0
 f: +49 (0)6021-58830-11
 anfrage@phenomenex.com

India

t: +91 (0)40-3012 2400
 f: +91 (0)40-3012 2411
 indiainfo@phenomenex.com

Ireland

t: +353 (0)1 247 5405
 f: +44 1625-501796
 eireinfo@phenomenex.com

Italy

t: +39 051 6327511
 f: +39 051 6327555
 italiainfo@phenomenex.com

Luxembourg

t: +31 (0)30-2418700
 f: +31 (0)30-2383749
 nlinfo@phenomenex.com

Mexico

t: 01-800-844-5226
 f: 001-310-328-7768
 tecnicomx@phenomenex.com

The Netherlands

t: +31 (0)30-2418700
 f: +31 (0)30-2383749
 nlinfo@phenomenex.com

New Zealand

t: +64 (0)9-4780951
 f: +64 (0)9-4780952
 nzinfo@phenomenex.com

Norway

t: +47 810 02 005
 f: +45 4810 6265
 nordicinfo@phenomenex.com

Portugal

t: +351 221 450 488
 f: +34 91-413-2290
 ptinfo@phenomenex.com

Spain

t: +34 91-413-8613
 f: +34 91-413-2290
 espinfo@phenomenex.com

Sweden

t: +46 (0)8 611 6950
 f: +45 4810 6265
 nordicinfo@phenomenex.com

Switzerland

t: +41 61 692 20 20
 f: +41 61 692 20 22
 swissinfo@phenomenex.com

United Kingdom

t: +44 (0)1625-501367
 f: +44 (0)1625-501796
 ukinfo@phenomenex.com

USA

t: +1 (310) 212-0555
 f: +1 (310) 328-7768
 info@phenomenex.com

All other countries Corporate Office USA

t: +1 (310) 212-0555
 f: +1 (310) 328-7768
 info@phenomenex.com

Terms and Conditions

Subject to Phenomenex Standard Terms & Conditions, which may be viewed at www.phenomenex.com/TermsAndConditions.

Trademarks

bioZen is a trademark of Phenomenex.

FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures.

© 2018 Phenomenex, Inc. All rights reserved.

