

HPLC Application Fast Analysis of Sugars in Fruit Juices - Moscato (White Wine) Rezex[™] RCM-Monosaccharide Ca+2 (8%), LC Column 100 x 7.8 mm, Ea Column: **Dimensions:** 100 x 7.8 mm ID ēΧ Order No: 00D-0130-K0 Elution Type: Isocratic Eluent A: Water Gradient Time (min) Step No. Pct A **Profile:** 0 100 1 Flow Rate: 0.4 mL/min Col. Temp.: 80 °C Products used in this applic **Detection:** Refractive Index (RI) @ (35 °C) 23056 **REZEX** nRIU EtOH 17500 3 15000 12500 10000 7500 2 5000

ANALYTES:

1 Moscato (White Wine)

2500

0

Ó

2

4

6

8

10

12

14 min

©2025 Phenomenex Inc. All rights reserved.

Phenomenex products are available worldwide.

www.phenomenex.com

For more information contact your Phenomenex Representative at support@phenomenex.com