
HPLC Application

Synthetic Cannabinoids by LCMS using Luna Omega 1.6u C18 50x2.1mm



©2025 Phenomenex Inc. All rights reserved.

Phenomenex products are available worldwide.

www.phenomenex.com

For more information contact your Phenomenex Representative at support@phenomenex.com



HPLC Application



Synthetic Cannabinoids by LCMS using Luna Omega 1.6u C18 50x2.1mm

ANALYTES:

- JWH-073 Butanoic acid Metabolite 1
- 2 JWH-073 3-hydroxybutyl metabolite (isomer to #1)
- 3 JWH-018 Pentanoic acid metabolite
- JWH-073 4-hydroxybutyl metabolite 4
- 5 JWH-018 4-hydroxypentyl metabolite (isomer to #3)
- 6 AM-2201 4-hydroxypentyl metabolite
- 7 JWH-018 5-hydroxypentyl metabolite (isomer to #3)
- 8 AM-694
- 9 AM-2201
- 10 JWH-073
- 11 JWH-018

©2025 Phenomenex Inc. All rights reserved.

Phenomenex products are available worldwide.

www.phenomenex.com

For more information contact your Phenomenex Representative at support@phenomenex.com