

EPA 8082 Aroclor 1260 by GC/ECD on ZB-CLPesticides-1 (Helium Gas)

Column: Zebtron ZB-CLPesticides-1 Column 30m x 0.25mm x 0.25um , Ea
Phase: Proprietary Pesticides Phase
Dimensions: 30 meters x 0.25 mm x 0.25 µm
Order No: 7HG-G028-11
Oven Profile: 120 °C to 200 @ 45 °C/min to 330 °C @ 15 °C/min for 2 min


Carrier Gas: Constant Flow Helium, 60 cm/sec
Injection: Splitless :1 1 µL @ 250°C
Detection: Electron Capture (ECD) (330°C)
Analyst Note: 1 µL Pulsed Splitless @ 30 psi for 20 sec.
Aroclor @ 1 ppm and surrogate and internal standard @ 100 ppb in hexane



Products used in this application:



22642



ANALYTES:

- 1 Tetrachloro-meta-xylene [TCMX] (Surr)
- 2 Decachlorobiphenyl (IS)

