


## Batch To Batch Reproducibility on bioZen Intact C4

**Column:** N/A  
**Dimensions:** N/A  
**Order No:** N/A  
**Elution Type:** N/A  
**Eluent A:** 0.1% TFA in Water  
**Eluent B:** 0.1% TFA in ACN

Gradient Profile:	Step No.	Time (min)	Pct A	Pct B
	1	0	90	10
	2	0	70	30
	3	0	50	50
	4	0	90	10

**Flow Rate:** 1.5 mL/min  
**Col. Temp.:** 40 °C  
**Detection:** UV-Vis Abs.-Variable Wave.(UV) @ 280 nm (ambient)



### ANALYTES:

- 1 Ligh Chain
- 2 Heavy Chain
- 3 Intact mAb

