


Batch To Batch Reproducibility on bioZen Intact C4

Column: N/A
Dimensions: N/A
Order No: N/A
Elution Type: N/A
Eluent A: 0.1% TFA in Water
Eluent B: 0.1% TFA in ACN

Gradient Profile:	Step No.	Time (min)	Pct A	Pct B
	1	0	90	10
	2	0	70	30
	3	0	50	50
	4	0	90	10

Flow Rate: 1.5 mL/min
Col. Temp.: 40 °C
Detection: UV-Vis Abs.-Variable Wave.(UV) @ 280 nm (ambient)



ANALYTES:

- 1 Ligh Chain
- 2 Heavy Chain
- 3 Intact mAb

