

Thioridazine HCl and Mesoridazine Besylate

Column: Ultracarb[™] 5 μ m ODS (30) 60 Å, LC Column 150 x 4.6 mm, Ea

Dimensions: 150 x 4.6 mm ID

Order No: 00F-0351-E0

Elution Type: Isocratic

Eluent A: Acetonitrile

Eluent B: Water with 0.01% TEA

Gradient Profile:	Step No.	Time (min)	Pct A	Pct B
	1	0	85	15

Flow Rate: 2.5 mL/min


Col. Temp.: ambient

Detection: UV-Vis Abs.-Variable Wave.(UV) @ 265 nm (ambient)

Analyst Note: SecurityGuard[™] Guard Cartridge System extends column lifetime.

- SecurityGuard Cartridges, C18 4 x 3.0mm, 10/Pk Part No.: AJ0-4287

- Holder Part No.: KJ0-4282



ANALYTES:

1 Mesoridazine besylate

2 Thioridazine



Products used in this application:

Ultracarb[™]

